

Glutamyl aminopeptidase

Catalog No: #AP80512



Package Size: #AP80512-1 50ug #AP80512-2 100ug #AP80512-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glutamyl aminopeptidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	EAP,Aminopeptidase A,AP-A,Differentiation antigen gp160,CD249
Accession No.	Uniprot ID:Q07075
Uniprot	Q07075
GeneID	2028;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ENPEP

Note: This product is for in vitro research use only