

Natriuretic peptides B

Catalog No: #AP80513



Package Size: #AP80513-1 50ug #AP80513-2 100ug #AP80513-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Natriuretic peptides B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Gamma-brain natriuretic peptide
Accession No.	Uniprot ID:P40753
Uniprot	P40753
GeneID	18158;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Nppb

Note: This product is for in vitro research use only