

# Annexin A4

Catalog No: #AP80514



Package Size: #AP80514-1 50ug #AP80514-2 100ug #AP80514-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Annexin A4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	36 kDa zymogen granule membrane-associated protein,ZAP36,Annexin IV,Annexin-4,Lipocortin IV,ANX4
Accession No.	Uniprot ID:P50994
Uniprot	P50994
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ANXA4

Note: This product is for in vitro research use only