

Annexin A6

Catalog No: #AP80516



Package Size: #AP80516-1 50ug #AP80516-2 100ug #AP80516-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Annexin A6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	67 kDa calelectrin,Annexin VI,Annexin-6,Calphobindin-II,CPB-II,Chromobindin-20,Lipocortin VI,P68,P70,Protein III,ANX6
Accession No.	Uniprot ID:P51901
Uniprot	P51901
GeneID	395481;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ANXA6

Note: This product is for in vitro research use only