

Aminoacyl tRNA synthase complex-interacting multifunctional protein 1

Catalog No: #AP80518

Package Size: #AP80518-1 50ug #AP80518-2 100ug #AP80518-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Aminoacyl tRNA synthase complex-interacting multifunctional protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Multisynthase complex auxiliary component p43,EMAP2,SCYE1
Accession No.	Uniprot ID:Q12904
Uniprot	Q12904
GeneID	9255;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AIMP1

Note: This product is for in vitro research use only