

Transferrin Monoclonal Antibody(10G3)

Catalog No: #44109

Package Size: #44109-1 50ul #44109-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	Transferrin Monoclonal Antibody(10G3)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	10G3
Isotype	IgG1
Applications	WB IHC
Species Reactivity	Hu
Specificity	Transferrin Mouse Monoclonal AntibodyB£B"10G3B£B© detects endogenous levels of Transferrin
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Target Name	Transferrin
Accession No.	Swiss-Prot:P02787Gene ID:7018
Uniprot	P02787
GeneID	7018;
Calculated MW	77kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western blotting: 1:2000-1:5000

Immunohistochemistry: 1:100-1:200


Images



Western blot analysis of Human Serum with Transferrin Mouse mAb diluted at 1:2,000.



Immunohistochemical analysis of paraffin-embedded human-liver using antibody diluted at 1:50.



IHC: 1:100-1:200

Note: This product is for in vitro research use only