

Alanine aminotransferase 2

Catalog No: #AP80521



Package Size: #AP80521-1 50ug #AP80521-2 100ug #AP80521-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Alanine aminotransferase 2 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | ALT2,Glutamate pyruvate transaminase 2,GPT 2,Glutamic--alanine transaminase 2,Glutamic--pyruvic transaminase 2,AAT2,ALT2 |
| Accession No. | Uniprot ID:Q8TD30 |
| Uniprot | Q8TD30 |
| GeneID | 84706; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GPT2

Note: This product is for in vitro research use only