

Adipocyte enhancer-binding protein 1

Catalog No: #AP80522



Package Size: #AP80522-1 50ug #AP80522-2 100ug #AP80522-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Adipocyte enhancer-binding protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	AE-binding protein 1,Aortic carboxypeptidase-like protein,Aclp
Accession No.	Uniprot ID:A2RUV9
Uniprot	A2RUV9
GeneID	305494;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Aebp1

Note: This product is for in vitro research use only