

Adenosine deaminase domain-containing protein 1

Catalog No: #AP80525



Package Size: #AP80525-1 50ug #AP80525-2 100ug #AP80525-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Adenosine deaminase domain-containing protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Testis nuclear RNA-binding protein,Tenr
Accession No.	Uniprot ID:Q5SUE7
Uniprot	Q5SUE7
GeneID	21744;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Adad1

Note: This product is for in vitro research use only