

Apelin receptor

Catalog No: #AP80526

Package Size: #AP80526-1 50ug #AP80526-2 100ug #AP80526-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Apelin receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Angiotensin receptor-like 1,B78,G-protein coupled receptor APJ,GPCR34,Agtr11,Apj
Accession No.	Uniprot ID:Q9JHG3
Uniprot	Q9JHG3
GeneID	83518;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ap1nr

Note: This product is for in vitro research use only