

Apolipoprotein A-I

Catalog No: #AP80527



Package Size: #AP80527-1 50ug #AP80527-2 100ug #AP80527-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Apolipoprotein A-I
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Apo-AI,ApoA-I,Apolipoprotein A1
Accession No.	Uniprot ID:P02648
Uniprot	P02648
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:APOA1

Note: This product is for in vitro research use only