

Apolipoprotein B-100

Catalog No: #AP80528



Package Size: #AP80528-1 50ug #AP80528-2 100ug #AP80528-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Apolipoprotein B-100
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Apo B-100,Aa1064,Ac1-060
Accession No.	Uniprot ID:Q7TMA5
Uniprot	Q7TMA5
GeneID	54225;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Apob

Note: This product is for in vitro research use only