

Transferrin Monoclonal Antibody(3A5)

Catalog No: #44110

Package Size: #44110-1 50ul #44110-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Transferrin Monoclonal Antibody(3A5)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	3A5
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu
Specificity	Transferrin Mouse Monoclonal AntibodyB£B™3A5B£B© detects endogenous levels of Transferrin
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Target Name	Transferrin
Accession No.	Swiss-Prot:P02787Gene ID:7018
Uniprot	P02787
GeneID	7018;
Calculated MW	77kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded rat-liver using antibody diluted at 1:50.

Note: This product is for in vitro research use only