

Single-stranded DNA cytosine deaminase

Catalog No: #AP80536



Package Size: #AP80536-1 50ug #AP80536-2 100ug #AP80536-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Single-stranded DNA cytosine deaminase |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | Activation-induced cytidine deaminase,AID,Cytidine aminohydrolase,Aid |
| Accession No. | Uniprot ID:Q9WVE0 |
| Uniprot | Q9WVE0 |
| GeneID | 11628; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Aicda

Note: This product is for in vitro research use only