

## Cytosol aminopeptidase

Catalog No: #AP80537



Package Size: #AP80537-1 50ug #AP80537-2 100ug #AP80537-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Cytosol aminopeptidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Leucine aminopeptidase 3,LAP-3,Leucyl aminopeptidase,Peptidase S,Proline aminopeptidase,Prolyl aminopeptidase,LAPEP,PEPS
Accession No.	Uniprot ID:P28838
Uniprot	P28838
GeneID	51056;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:LAP3

Note: This product is for in vitro research use only