

Angiotensin-converting enzyme 2

Catalog No: #AP80538



Package Size: #AP80538-1 50ug #AP80538-2 100ug #AP80538-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Angiotensin-converting enzyme 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	425O18-1,Angiotensin-like protein 4,Fasting-induced adipose factor,Hepatic fibrinogen/angiotensin-related protein,HFARP,Secreted protein Bk89,Farp,Fiaf,Ng27
Accession No.	Uniprot ID:Q9Z1P8
Uniprot	Q9Z1P8
GeneID	57875;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Angptl4

Note: This product is for in vitro research use only