

Apoptotic protease-activating factor 1

Catalog No: #AP80539



Package Size: #AP80539-1 50ug #AP80539-2 100ug #AP80539-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Apoptotic protease-activating factor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	APAF-1
Accession No.	Uniprot ID:Q9GL23
Uniprot	Q9GL23
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:APAF1

Note: This product is for in vitro research use only