

# Ankyrin-2

Catalog No: #AP80542



Package Size: #AP80542-1 50ug #AP80542-2 100ug #AP80542-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Ankyrin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ANK-2,Ankyrin-B,Brain ankyrin
Accession No.	Uniprot ID:Q8C8R3
Uniprot	Q8C8R3
GeneID	109676;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Ank2

Note: This product is for in vitro research use only