

Angiotensin-converting enzyme 3

Catalog No: #AP80544



Package Size: #AP80544-1 50ug #AP80544-2 100ug #AP80544-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Angiotensin-converting enzyme 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Angiotensin-converting enzyme 3
Accession No.	Uniprot ID:Q9R182
Uniprot	Q9R182
GeneID	30924;
SDS-PAGE MW	22.88Kda
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:93-300aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Angptl3

Note: This product is for in vitro research use only