

C-type natriuretic peptide

Catalog No: #AP80549



Package Size: #AP80549-1 50ug #AP80549-2 100ug #AP80549-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	C-type natriuretic peptide
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	CNP
Accession No.	Uniprot ID:P18104
Uniprot	P18104
GeneID	493772;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: NPPC

Note: This product is for in vitro research use only