

Galectin-3 Monoclonal Antibody(6B8)

Catalog No: #44113

Package Size: #44113-1 50ul #44113-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Galectin-3 Monoclonal Antibody(6B8)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	6B8
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu
Specificity	Galectin-3 Mouse Monoclonal Antibody(6B8) detects endogenous levels of Galectin-3
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Target Name	Galectin-3
Accession No.	Swiss-Prot:P17931Gene ID:3958
Uniprot	P17931
GeneID	3958;
Calculated MW	26kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded human-colon using antibody diluted at 1:50.

Note: This product is for in vitro research use only