

Aminopeptidase N

Catalog No: #AP80556



Package Size: #AP80556-1 50ug #AP80556-2 100ug #AP80556-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Aminopeptidase N
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	AP-N,Alanyl aminopeptidase,Aminopeptidase M,AP-M,Kidney Zn peptidase,KZP,Microsomal aminopeptidase,CD13
Accession No.	Uniprot ID:P15684
Uniprot	P15684
GeneID	81641;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Anpep

Note: This product is for in vitro research use only