

Transcription factor AP-2-beta

Catalog No: #AP80558



Package Size: #AP80558-1 50ug #AP80558-2 100ug #AP80558-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Transcription factor AP-2-beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	AP2-beta,Activating enhancer-binding protein 2-beta
Accession No.	Uniprot ID:Q92481
Uniprot	Q92481
GeneID	7021;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TFAP2B

Note: This product is for in vitro research use only