

Galectin-3 Monoclonal Antibody(8D7)

Catalog No: #44114

Package Size: #44114-1 50ul #44114-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	Galectin-3 Monoclonal Antibody(8D7)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	8D7
Isotype	IgG1
Applications	WB IHC
Species Reactivity	Hu Ms
Specificity	Galectin-3 Mouse Monoclonal Antibody(8D7) detects endogenous levels of Galectin-3
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Target Name	Galectin-3
Accession No.	Swiss-Prot:P17931Gene ID:3958
Uniprot	P17931
GeneID	3958;
Calculated MW	26kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western blotting: 1:2000-1:5000

Immunohistochemistry: 1:100-1:200


Images



Western blot analysis of 1)MCF7, 2) 3T3 with Galectin-3 Mouse mAb diluted at 1:2,000.



Immunohistochemical analysis of paraffin-embedded human-colon using antibody diluted at 1:50



IHC: 1:100-1:200

Note: This product is for in vitro research use only