

## Aminopeptidase N

Catalog No: #AP80561



Package Size: #AP80561-1 50ug #AP80561-2 100ug #AP80561-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Aminopeptidase N
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	AP-N,cAPN,Alanyl aminopeptidase,Aminopeptidase M,AP-M,Microsomal aminopeptidase,CD13
Accession No.	Uniprot ID:P79143
Uniprot	P79143
GeneID	403913;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ANPEP

Note: This product is for in vitro research use only