

# Pancreatic alpha-amylase

Catalog No: #AP80563



Package Size: #AP80563-1 50ug #AP80563-2 100ug #AP80563-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Pancreatic alpha-amylase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PA, 1,4-alpha-D-glucan glucanohydrolase, Amy2a
Accession No.	Uniprot ID:P00688
Uniprot	P00688
GeneID	100043684;100043686;100043688;109959;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Amy2

Note: This product is for in vitro research use only