

## Membrane primary amine oxidase

Catalog No: #AP80564



Package Size: #AP80564-1 50ug #AP80564-2 100ug #AP80564-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Membrane primary amine oxidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Copper amine oxidase,HPAO,Semicarbazide-sensitive amine oxidase,SSAO,Vascular adhesion protein 1,VAP-1,VAP1
Accession No.	Uniprot ID:Q16853
Uniprot	Q16853
GeneID	8639;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:AOC3

Note: This product is for in vitro research use only