

Acetylserotonin O-methyltransferase

Catalog No: #AP80566



Package Size: #AP80566-1 50ug #AP80566-2 100ug #AP80566-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Acetylserotonin O-methyltransferase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Hydroxyindole O-methyltransferase,Hiomt
Accession No.	Uniprot ID:B3GSH5
Uniprot	B3GSH5
GeneID	246281;
Tag Info	His
Formulation	Storage Buffer:20mM Tris,0.5mM EDTA,0.5M Arginine,20%glycerol (PH8.0)
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)

gene_full_name:Asmt

Note: This product is for in vitro research use only