

Alanine aminotransferase 1

Catalog No: #AP80568



Package Size: #AP80568-1 50ug #AP80568-2 100ug #AP80568-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Alanine aminotransferase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	ALT1,Glutamate pyruvate transaminase 1,GPT 1,Glutamic--alanine transaminase 1,Glutamic--pyruvic transaminase 1,GPT1
Accession No.	Uniprot ID:A4IFH5
Uniprot	A4IFH5
GeneID	539188;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GPT

Note: This product is for in vitro research use only