

Aldehyde dehydrogenase, mitochondrial

Catalog No: #AP80571



Package Size: #AP80571-1 50ug #AP80571-2 100ug #AP80571-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Aldehyde dehydrogenase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	ALDH class 2,ALDH-E2
Accession No.	Uniprot ID:Q2XQV4
Uniprot	Q2XQV4
GeneID	733685;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ALDH2

Note: This product is for in vitro research use only