

E3 ubiquitin-protein ligase AMFR

Catalog No: #AP80573



Package Size: #AP80573-1 50ug #AP80573-2 100ug #AP80573-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	E3 ubiquitin-protein ligase AMFR
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Autocrine motility factor receptor, isoform 2,AMF receptor, isoform 2,RING finger protein 45,gp78,RNF45
Accession No.	Uniprot ID:Q9UKV5
Uniprot	Q9UKV5
GeneID	267;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AMFR

Note: This product is for in vitro research use only