

Anti-Muellerian hormone type-2 receptor

Catalog No: #AP80577



Package Size: #AP80577-1 50ug #AP80577-2 100ug #AP80577-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Anti-Muellerian hormone type-2 receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Anti-Muellerian hormone type II receptor,AMH type II receptor,MIS type II receptor,MISRII
Accession No.	Uniprot ID:Q8K592
Uniprot	Q8K592
GeneID	110542;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Amhr2

Note: This product is for in vitro research use only