

## Transferrin Monoclonal Antibody(5E5)

Catalog No: #44116

Package Size: #44116-1 50ul #44116-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Transferrin Monoclonal Antibody(5E5)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	5.00E+05
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu
Specificity	Transferrin Mouse Monoclonal AntibodyB£B`5E5£B© detects endogenous levels of Transferrin
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Target Name	Transferrin
Accession No.	Swiss-Prot:P02787Gene ID:7018
Uniprot	P02787
GeneID	7018;
Calculated MW	77kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Immunohistochemistry: 1:100-1:200

## Images



Immunohistochemical analysis of paraffin-embedded rat-liver using antibody diluted at 1:50.

---

Note: This product is for in vitro research use only