

Adenomatous polyposis coli protein

Catalog No: #AP80580



Package Size: #AP80580-1 50ug #AP80580-2 100ug #AP80580-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Adenomatous polyposis coli protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Protein APC, Deleted in polyposis 2.5, DP2.5
Accession No.	Uniprot ID:P25054
Uniprot	P25054
GeneID	324;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: APC

Note: This product is for in vitro research use only