

Aquaporin-4

Catalog No: #AP80582



Package Size: #AP80582-1 50ug #AP80582-2 100ug #AP80582-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Aquaporin-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	AQP-4,Mercurial-insensitive water channel,MIWC,WCH4
Accession No.	Uniprot ID:P47863
Uniprot	P47863
GeneID	25293;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Aqp4

Note: This product is for in vitro research use only