

Antithrombin-III

Catalog No: #AP80583



Package Size: #AP80583-1 50ug #AP80583-2 100ug #AP80583-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Antithrombin-III
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ATIII,Serpin C1,At3
Accession No.	Uniprot ID:P32261
Uniprot	P32261
GeneID	11905;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Serpinc1

Note: This product is for in vitro research use only