

## Annexin A3

Catalog No: #AP80584



Package Size: #AP80584-1 50ug #AP80584-2 100ug #AP80584-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Annexin A3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	35-alpha calcimedlin,Annexin III,Annexin-3,Inositol 1,2-cyclic phosphate 2-phosphohydrolase,Lipocortin III,Placental anticoagulant protein III,PAP-III,ANX3
Accession No.	Uniprot ID:P12429
Uniprot	P12429
GeneID	306;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ANXA3

Note: This product is for in vitro research use only