

Annexin A5

Catalog No: #AP80586



Package Size: #AP80586-1 50ug #AP80586-2 100ug #AP80586-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Annexin A5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Anchorin CII,Annexin V,Annexin-5,Calphobindin I,CBP-I,Endonexin II,Lipocortin V,Placental anticoagulant protein 4,PP4,Placental anticoagulant protein I,PAP-I,Thromboplastin inhibitor,Vascular anticoagulant-alpha,VAC-alpha,Anx5
Accession No.	Uniprot ID:P14668
Uniprot	P14668
GeneID	25673;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Anxa5

Note: This product is for in vitro research use only