

Rho guanine nucleotide exchange factor 1

Catalog No: #AP80589



Package Size: #AP80589-1 50ug #AP80589-2 100ug #AP80589-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rho guanine nucleotide exchange factor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	115 kDa guanine nucleotide exchange factor,p115-RhoGEF,Sub1.5
Accession No.	Uniprot ID:Q92888
Uniprot	Q92888
GeneID	9138;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ARHGEF1

Note: This product is for in vitro research use only