

Annexin A1

Catalog No: #AP80590



Package Size: #AP80590-1 50ug #AP80590-2 100ug #AP80590-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Annexin A1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Annexin I,Annexin-1,Calpactin II,Calpactin-2,Chromobindin-9,Lipocortin I,Phospholipase A2 inhibitory protein,p35,Anx1,Lpc-1,Lpc1
Accession No.	Uniprot ID:P10107
Uniprot	P10107
GeneID	16952;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Anxa1

Note: This product is for in vitro research use only