

Androgen receptor

Catalog No: #AP80594



Package Size: #AP80594-1 50ug #AP80594-2 100ug #AP80594-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Androgen receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Dihydrotestosterone receptor,Nuclear receptor subfamily 3 group C member 4,NR3C4
Accession No.	Uniprot ID:Q9TT90
Uniprot	Q9TT90
GeneID	403588;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AR

Note: This product is for in vitro research use only