

Natriuretic peptides A

Catalog No: #AP80595



Package Size: #AP80595-1 50ug #AP80595-2 100ug #AP80595-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Natriuretic peptides A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CDD-ANF,Cardiodilatin,CDD,Cardiodilatin-related peptide,CDP,Prepronatriodilatin,ANP,PND
Accession No.	Uniprot ID:P01160
Uniprot	P01160
GeneID	4878;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NPPA

Note: This product is for in vitro research use only