

Fructose-bisphosphate aldolase C

Catalog No: #AP80597



Package Size: #AP80597-1 50ug #AP80597-2 100ug #AP80597-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fructose-bisphosphate aldolase C
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Aldolase 3,Brain-type aldolase,Scrapie-responsive protein 2,Zebrin II,Aldo3,Scrg2
Accession No.	Uniprot ID:P05063
Uniprot	P05063
GeneID	11676;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Aldoc

Note: This product is for in vitro research use only