

Type-1 angiotensin II receptor

Catalog No: #AP80604



Package Size: #AP80604-1 50ug #AP80604-2 100ug #AP80604-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Type-1 angiotensin II receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	AT1AR,AT1BR,Angiotensin II type-1 receptor,AT1,AGTR1A,AGTR1B,AT2R1,AT2R1B
Accession No.	Uniprot ID:P30556
Uniprot	P30556
GeneID	185;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AGTR1

Note: This product is for in vitro research use only