

Autoimmune regulator

Catalog No: #AP80605



Package Size: #AP80605-1 50ug #AP80605-2 100ug #AP80605-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Autoimmune regulator
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Autoimmune polyendocrinopathy candidiasis ectodermal dystrophy protein homolog,APECED protein homolog
Accession No.	Uniprot ID:Q9Z0E3
Uniprot	Q9Z0E3
GeneID	11634;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Aire

Note: This product is for in vitro research use only