

Transcription factor AP-2-beta

Catalog No: #AP80608



Package Size: #AP80608-1 50ug #AP80608-2 100ug #AP80608-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Transcription factor AP-2-beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	AP2-beta,Activating enhancer-binding protein 2-beta,Tcfap2b
Accession No.	Uniprot ID:Q61313
Uniprot	Q61313
GeneID	21419;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tfap2b

Note: This product is for in vitro research use only