

# Transcription factor AP-2-beta

Catalog No: #AP80612



Package Size: #AP80612-1 50ug #AP80612-2 100ug #AP80612-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Transcription factor AP-2-beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	AP2-beta,Activating enhancer-binding protein 2-beta
Accession No.	Uniprot ID:Q76HI7
Uniprot	Q76HI7
GeneID	403463;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:TFAP2B

Note: This product is for in vitro research use only