

Apelin receptor

Catalog No: #AP80613



Package Size: #AP80613-1 50ug #AP80613-2 100ug #AP80613-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Apelin receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Angiotensin receptor-like 1,G-protein coupled receptor APJ,MSR,Agtr1,Apj
Accession No.	Uniprot ID:Q9WV08
Uniprot	Q9WV08
GeneID	23796;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Aplnr

Note: This product is for in vitro research use only