

Angiopoietin-4

Catalog No: #AP80615



Package Size: #AP80615-1 50ug #AP80615-2 100ug #AP80615-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Angiopoietin-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ANG-4,Angiopoietin-3,ANG-3,ANG3,ANG4
Accession No.	Uniprot ID:Q9Y264
Uniprot	Q9Y264
GeneID	51378;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ANGPT4

Note: This product is for in vitro research use only