

Angiotensin-converting enzyme 2

Catalog No: #AP80627

Package Size: #AP80627-1 50ug #AP80627-2 100ug #AP80627-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Angiotensin-converting enzyme 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Angiotensin-converting enzyme 2
Accession No.	Uniprot ID:Q2KJ51
Uniprot	Q2KJ51
GeneID	509963;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ANGPTL4

Note: This product is for in vitro research use only